Nalone Limited

Web: www.nalone.co.uk
E-mail: info@nalone.co.uk

Yany **USER MANUAL**

Yany

Language

➤ English	-----	02-03
➤ Français	-----	04-05
➤ Italiano	-----	06-07
➤ Português	-----	08-09
➤ Русский	-----	10-11
➤ Español	-----	12-13
➤ Svenska	-----	14-15
➤ Deutsch	-----	16-18

We thank you for purchasing our products.
Please read this manual carefully before use!


➤ English

Specifications

- Ball Material:Body safe metal
- Ball Retainer:Body safe silicone
- Gold ball weight:30g
- Blue ball weight:40g
- Ball diameter:31mm
- Overall dimensions:155mm(L) x32mm(D)

How to use

1. Clean before use
2. Use water-based lubricant before use
3. Step 1: Empty your bladder before inserting the balls into the vagina.
Step 2: One ball should be inserted at a time. Inserting them is similar to inserting vaginal tampons. In case you have trouble inserting them, lift one leg up or insert them while lying down. Keep the silicone retainer cord outside your body for ball removal.
Step 3: To hold the balls in the vagina, squeeze first the muscles of your leg and then your PC muscles. See diagram below:
Step 4: Kegel balls should be held inside your vagina for at least 15 minutes in a day so that your PC muscles are strengthened. For a better workout you can hold them inside for Longer.


- After finishing your exerciser lie down, relax and pull on the silicone retainer tail.
- After removing separate the Silicone retainer and metal balls and carefully clean. Please ensure you wash/disinfect you massager with antibacterial soap and warm water.

Training Tips

Always begin using the lighter ball and exercise for 15 minutes. You can add 5 minutes each day until you are exercising for 30 minutes per session.

Suggested training plan

- Week 1 - 1 30g golden ball
- Week 2 - 1 40g blue ball
- Week 3 & 4 - Link 2 30g golden balls
- Week 5 & 6 - Link 1 30g golden ball and 1 40g blue ball
- Week 7 & 8 - Link 2 40g golden balls

Warnings

- ◆ Please do not use during Pregnancy or Menstrual cycle
- ◆ You may feel some discomfort the first day after use, as your body gets use to the fiscal exercise.
- ◆ Should this discomfort continue or you experience any other side effects please stop using and consult your doctor.

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   
NALONE ALL RIGHTS RESERVED

◆ What is Covered

This warranty covers any defects in the Nalone product purchased provided the product has been used and maintained in accordance with terms of the Nalone instruction manual.

◆ How Long Coverage Lasts

This warranty runs for one year from the date of purchase of your Nalone product. Please retain your receipt as proof of purchase.

◆ What is not Covered

This warranty does not cover accidental or wilful damage or any defects arising from the product not being cared for as per the instruction manual. Please take care to carefully clean and store the product.

◆ What Nalone Will Do

Nalone will repair any product that proves to be defective in materials or workmanship. In the event repair is not possible, Nalone will either replace, or refund the full purchase price, whichever you prefer.

◆ How To Get Service

Contact Nalone on info@nalone.co.uk

User Information

User name:	_____	Tel:	_____
User Address:	_____		
Product name:	_____		
Purchase date:	_____		

About Us

Nalone is a world leader in the design and production of adult lifestyle products. The company's philosophy has been to devise and develop a premium range of intimate massagers. The in-house design team are constantly looking to innovate and strive towards creating the ultimate user experience.

All products are presented in exquisite packaging and where necessary include a USB charger, storage bag and instruction for use and product care.

All products carry Nalone's standard 1 year warranty.

www.nalone.co.uk

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Nalone Ltd.